


Aspartate Aminotransferase (AST)

Test principle: UV test



Method according to the International Federation of Clinical Chemistry (IFCC), without pyridoxal-5'-phosphate. AST in the sample catalyzes the transfer of an amino group between L-aspartate and α -ketoglutarate to form oxaloacetate and L-glutamate. The oxaloacetate then reacts with NADH in the presence of malate dehydrogenase (MDH) to form NAD^+ .

The rate of the NADH oxidation is directly proportional to the catalytic AST activity. It is determined by measuring the decrease in absorbance at 340 nm.

Proposed reagent composition

approximately 2+1 formulation

Reagent 1

Composition	Concentration	Catalog Number
Buffer (TRIS, pH 7.8)	280 mmol/l	10 153 265 001
L-Aspartate	800 mmol/l	
Malate dehydrogenase (MDH)	>2.1 kU/l	10 267 155 103
Lactate dehydrogenase (LDH)	>3.0 kU/l	10 679 666 103 or 12 235 650 103
Detergent, preservative, stabilizer, such as Sodium azide		
Triton X-100	0.01 %	10 743 119 103
Albumin	0.25 %	10 738 328 103

Reagent 2

Composition	Concentration	Catalog Number
NADH	>1.7 mmol/l	10 004 642 103
α -Ketoglutarate, di-Na	94 mmol/l	10 040 584 103
Preservative, such as Sodium azide		

Products are for further processing only.

All brands or product names are trademarks of their respective holders.

custombiotech.roche.com

Your Roche Custom Biotech Customer Service

Europe, Middle East, Africa, Latin America

Phone +49 621 759 8580

Fax +49 621 759 8610

mannheim.custombiotech@roche.com

Japan Phone +81 3 5443 5285

Fax +81 3 5443 7934

japan.custombiotech@roche.com

Asia Pacific Phone +65 6371 6638

Fax +65 6371 6601

apac.custombiotech@roche.com

United States

Phone +1 800 428 5433, ext. 14649 (toll-free)

Fax +1 317 521 4065

custombiotech.ussales@roche.com

Canada Phone +1 450 686 7050

Fax +1 450 686 7012

custombiotech.can@roche.com

Published by

Roche Diagnostics GmbH

Sandhofer Straße 116

68305 Mannheim

Germany

© 2011 Roche Diagnostics.

All rights reserved.

05837723990 © 0311