



Residual DNA *E. coli* Kit

High sensitivity for optimal product safety

The removal of host cell impurities is a critical step in the production of biopharmaceutical products. The Roche Residual DNA *E. coli* Kit is designed and validated for detection of residual *E. coli* DNA in biopharmaceutical products (protein solution, antibodies, etc.) as required by regulatory authorities.

Based on proven real-time qPCR technology, this kit makes detection of residual DNA from *E. coli* bacteria fast and reliable*.

It was developed to meet sensitivity requirements defined by WHO and FDA (10 ng *E. coli* DNA per therapeutic dose).

Achieve high sensitivity and specificity through improved tests*

- Linearity: 10 µg/ml to 5 pg/ml
- Lower quantitation limit: 5 pg/ml (50 fg/reaction)
- Lower detection limit: 1 pg/ml (10 fg/reaction)
- No cross-reactivity with unrelated DNA

Rely on a robust test procedure for a wide range of sample types

- Suitable for different matrices and a broad range of sample types from both in-process samples to bulk drug substances
- High lot-to-lot consistency of kit reagents, leading to reproducible results


Save time with proven real-time qPCR technology

- Time to result <5 hours*


Characteristics

Residual DNA *E. coli* Kit

Sample prep	Manual protocol optimized for excellent DNA recovery (min 80%).
Control concept	DNA stock solution for preparation of extraction controls and internal controls for PCR performance included in kit.
All reagents included	Supplied and ready-to-use reagents jump start your assay.
High quality	Roche is certified according to ISO 13485. Change notification available upon request.



01
 Typical analysis results obtained with Standard A (100 000 pg/ml) to G (1 pg/ml). Standard curve is calculated to Standard F (5 pg/ml LOQ) and should be linear. PCR efficiency should be 1.85 or higher.



02
 Typical results of a DNA spiked sample, a workflow positive and workflow negative controls. DNA recovery = minimum of 80% required.

Ordering information

Products

Catalog Number

Residual DNA *E. coli* Kit

07 728 735 001

QC Sample Preparation Kit

08 146 829 001

Additional Quality Control or In-Process Manufacturing Kits

MycoTOOL Mycoplasma Real-Time PCR Kit (160 PCR reactions)

06 495 605 001

MycoTOOL Mycoplasma Detection Amplification Kit

05 184 240 001

Residual DNA CHO Kit

07 427 689 001

Instrument

LightCycler® 480 Instrument II

05 015 278 001

Regulatory Disclaimer

For use in quality control/manufacturing process only.

The LightCycler® 480 Instrument is for life science research only. Not for use in diagnostic procedures.

Trademarks

LIGHTCYCLER, HIGH PURE, MYCOTOOL are trademarks of Roche. All other product names and trademarks are the property of their respective owners.

© 2023

All rights reserved.

Published by

Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany

[custombiotech.roche.com](https://www.custombiotech.roche.com)

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America

mannheim.custombiotech@roche.com

United States

custombiotech.ussales@roche.com

Canada

custombiotech.can@roche.com

Asia Pacific

apac.custombiotech@roche.com