


## Formate Assay for Cedex Bio and Cedex Bio HT Analyzers


*High performance and reproducibility with accurate data*

The Formate Assay developed for the Cedex Bio and Cedex Bio HT Analyzers provides fast and accurate quantitative measurements in cell culture and fermentation media. Formate is a metabolic intermediate relevant in the fermentation of bacteria and yeast. The Assay principle is based on the oxidation of formate to bicarbonate by NAD in the presence of formate dehydrogenase (FDH). The increase of NADH is directly proportional to the formate concentration and measured photometrically at 340 nm.



### Trust results to make high confidence decisions

- Accurate results over a wide linear measuring range
- Extended range with automated on-board dilution



**Figure 1:** Formate was spiked into RPMI cell culture medium and measured using the Cedex Bio HT Analyzer. (Verification data of Roche R&D)

### Benefit from a wide and sensitive measuring range

#### Formate Bio and Formate Bio HT

**Range:** 0.2 - 20 mmol/L, 9 - 900 mg/L  
up to 2 mol/L, 90 g/L with automated on-board dilution

### Conserve sample with low volume requirement

- Sample volume of  $\leq 5 \mu\text{l}$  per test

## Rely on precise results

	Level 1	Level 2	Level 3
<b>Mean</b>	1.0 mmol/L	8.0 mmol/L	15.0 mmol/L
<b>CV in-run</b>	0.8 %	0.3 %	0.4 %
<b>CV inter-run</b>	1.9 %	1.9 %	2.1 %

Precision was determined in samples of three concentration levels. Coefficients of variation (CV) were calculated for in-run precision (n=21) and inter-run precision (on 10 days). Representative performance data from Cedex Bio HT Analyzers are shown. Results obtained in individual laboratories may differ.

## Save time with improved workflow efficiency

- Ready-to-use reagents
- Long on-board and calibration curve stability
- Barcoded reagents
- Pre- and post-dilution capability

	Formate Bio	Formate Bio HT
On-board stability	21 days	140 days
Calibration interval	84 days	140 days

## Take control of your bioprocess with the Cedex portfolio

- Broad and expandable assay menu
- 2 instrument platforms with **identical technology, menu and reagents** to fit throughput and automation needs
- Solutions for development and manufacturing environments

## Ordering information

Product	Pack size	Catalog number
Formate Bio	4 x 50 tests	07 705 590 001
Formate Bio HT	200 tests	07 705 603 001
Calibrator D Bio	6 x 1 ml	07 368 321 001
Control D Level 1 Bio	6 x 1 ml	07 368 178 001
Control D Level 2 Bio	6 x 1 ml	07 368 186 001
Control D Level 3 Bio	6 x 1 ml	07 368 194 001

### Regulatory disclaimer

For use in quality control/manufacturing process only.

### Trademarks

CEDEX is a trademark of Roche.

[custombiotech.roche.com/cedex](http://custombiotech.roche.com/cedex)

## Please contact your local CustomBiotech representative

### Europe, Middle East, Africa, Latin America

Phone +49 621 759 8580  
Fax +49 621 759 6385  
[mannheim.custombiotech@roche.com](mailto:mannheim.custombiotech@roche.com)

### United States

Phone +1 800 428 5433, ext. 14649 (toll-free)  
Fax +1 317 521 4065  
[custombiotech.ussales@roche.com](mailto:custombiotech.ussales@roche.com)

### Canada

Phone +1 450 686 7050  
Fax +1 450 686 7012  
[custombiotech.can@roche.com](mailto:custombiotech.can@roche.com)

### Japan

Phone +81 3 6634 1046  
Fax +81 3 5479 0585  
[japan.custombiotech@roche.com](mailto:japan.custombiotech@roche.com)

### Asia Pacific

Phone +65 6371 6638  
Fax +65 6371 6601  
[apac.custombiotech@roche.com](mailto:apac.custombiotech@roche.com)

### Published by

Roche Diagnostics GmbH  
Sandhofer Straße 116  
68305 Mannheim  
Germany

© 2017 Roche Diagnostics GmbH  
All rights reserved.