



Residual DNA CHO Kit

High sensitivity for reliable detection

The removal of host cell impurities is a critical step in the production of biopharmaceutical products. The Residual DNA CHO kit is designed and validated for detection of residual CHO DNA in biopharmaceutical products (protein solution, antibodies, etc.) as required by regulatory authorities.

Based on proven real-time qPCR technology, this kit makes detection of residual DNA from Chinese Hamster Ovary (CHO) cell lines fast and reliable*. It was developed to meet sensitivity requirements defined by WHO and U.S. FDA (10 ng CHO DNA per therapeutic dose).

Achieve high sensitivity and specificity*

- Linearity: 4 pg/ml to 400 ng/ml
- Lower quantitation limit: 4 pg/ml (20 fg/reaction)
- Lower detection limit: 0.4 pg/ml (2 fg/reaction)
- No cross-reactivity with unrelated DNA

Rely on a robust test procedure for a wide range of sample types

- Suitable for different matrices and a broad range of sample types from both in-process samples to bulk drug substances
- High lot-to-lot consistency of kit reagents, leading to reproducible results



Save time with proven real-time qPCR technology

- Time to result <5 hours*

Characteristics


Residual DNA CHO Kit

Sample preparation	Manual protocol optimized for excellent DNA recovery (min 80%).
Control concept	DNA stock solution for preparation of extraction controls and internal controls for PCR performance included in kit.
All reagents included	Supplied and ready-to-use reagents jump start your assay.
High quality	Roche is certified according to ISO 13485. Change notification available upon request.



Error: 0,00522 | Efficiency: 1,948 | Skope: -3,454 | Yintercept: 36,20 | Link: 0,000

01
Typical analysis result with Standard A (400,000 pg/ml) to G (0.4 pg/ml). The standard curves should be linear and the PCR efficiency at least 1.85 or higher.



02
Typical analysis result for one sample, one sample spiked with DNA, one workflow positive control, and one workflow negative control. DNA recovery = 90%

Ordering information

Products

Catalog Number

Residual DNA CHO Kit

07 427 689 001

QC Sample Preparation Kit

08 146 829 001

Additional Quality Control Kits

MycoTOOL Mycoplasma Real-Time PCR Kit (160 PCR reactions)

06 495 605 001

MycoTOOL Mycoplasma Detection Amplification Kit

05 184 240 001

Residual DNA *E.coli* Kit

07 728 735 001

Instrument

LightCycler® 480 Instrument II

05 015 278 001

The Residual DNA CHO Kit can also be used with other real-time PCR platforms (e.g., QuantStudio™ 6 Flex PCR System from Applied Biosystems®).

Regulatory Disclaimer

For use in quality control/manufacturing process only.

The LightCycler® 480 Instrument is for life science research only. Not for use in diagnostic procedures.

Trademarks

LIGHTCYCLER, HIGH PURE, MYCOTOOL are trademarks of Roche. All other product names and trademarks are the property of their respective owners.

© 2025

All rights reserved.

Published by

Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany

[custombiotech.roche.com](https://www.custombiotech.roche.com)

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America

mannheim.custombiotech@roche.com

United States

custombiotech.ussales@roche.com

Canada

custombiotech.can@roche.com

Asia Pacific

apac.custombiotech@roche.com