

## Ethanol Assay for Cedex Bio and Cedex Bio HT Analyzers

### *Reliable determination with high specificity*

The Ethanol assay developed for the Cedex Bio and Cedex Bio HT Analyzers provides fast and accurate quantitative determination of ethyl alcohol (ethanol) in samples of fermentation media and other aqueous solutions.

#### Assay principle

The assay principle is based on the enzymatic oxidation of ethanol with NAD by a specific alcohol dehydrogenase (ADH). The resulting increase of NADH is measured photometrically at 340 nm.



#### Benefit from a wide measuring range


	Ethanol Bio and Ethanol Bio HT
<b>Range</b>	0.5 to 10 g/L, 11 to 220 mmol/L up to 200 g/L, 4.4 mol/L with automated sample predilution capability

#### Rely on high specificity

The assay uses an alcohol dehydrogenase enzyme variant with high specificity for ethanol. There is no influence of methanol. A low cross-reactivity was observed with propanol, butanol and pentanol, however, these signals are much lower compared to those of ethanol in similar concentration.

The following substances were tested for interference in the Ethanol assay using a high concentration of 20 g/L:

Substance (20 g/L)	Recovery
Methanol	negative (< 0.05 %)
<b>Ethanol</b>	<b>100 %</b>
1-Propanol	<b>15.1 %</b>
2-Propanol	negative (< 0.05 %)
1-Butanol	<b>4.3 %</b>
2-Methyl-1-propanol (Isobutanol)	negative (< 0.05 %)
1-Pentanol	<b>1.5 %</b>
Ethylene glycol	negative (< 0.05 %)
Glycerol	negative (< 0.05 %)
2-Phenoxyethanol (Phenylglycol)	negative (< 0.05 %)
Acetaldehyde	negative (< 0.05 %)



## Rely on precise results

	Level 1	Level 2	Level 3
Ethanol conc.	2.3 g/L	4.6 g/L	6.9 g/L
CV in-run	0.9 %	1.1 %	0.6 %
CV inter-run	1.2 %	1.1 %	0.9 %

Representative performance data from Cedex Bio HT Analyzers are shown. Results obtained in individual laboratories may differ. Coefficients of variation (CV) were calculated for in-run precision (n=21) and inter-run precision (on 10 days).

## Save time with improved workflow efficiency

- Fully automated assay
- No sample filtration or pretreatment required
- Ready-to-use reagents
- Barcoded reagents
- Long on-board and calibration curve stability:

	Ethanol Bio	Ethanol Bio HT
On-board stability	28 days	56 days
Calibration interval	56 days	56 days

## Take control of your bioprocess with the Cedex portfolio

- Broad and expandable assay menu
- Two instruments with identical technology, menu and reagents to meet your throughput and automation needs
- Appropriate solutions for development labs and qualified manufacturing environments

[custombiotech.roche.com/cedex](http://custombiotech.roche.com/cedex)

## Please contact your local CustomBiotech representative

### Europe, Middle East, Africa, Latin America

Phone +49 621 759 8580  
Fax +49 621 759 6385  
[mannheim.custombiotech@roche.com](mailto:mannheim.custombiotech@roche.com)

### United States

Phone +1 800 428 5433, ext. 14649 (toll-free)  
Fax +1 317 521 4065  
[custombiotech.ussales@roche.com](mailto:custombiotech.ussales@roche.com)

### Canada

Phone +1 450 686 7050  
Fax +1 450 686 7012  
[custombiotech.can@roche.com](mailto:custombiotech.can@roche.com)

### Japan

Phone +81 3 6634 1046  
Fax +81 3 5479 0585  
[japan.custombiotech@roche.com](mailto:japan.custombiotech@roche.com)

### Asia Pacific

Phone +65 6371 6638  
Fax +65 6371 6601  
[apac.custombiotech@roche.com](mailto:apac.custombiotech@roche.com)

## Ordering information

Product	Pack size	Catalog number
Ethanol Bio	4 x 50 tests	08 055 645 001
Ethanol Bio HT	200 tests	08 055 661 001
Calibrator E Bio	6 x 1 ml	08 083 703 001
Control E Level 1 Bio	6 x 1 ml	08 083 797 001
Control E Level 2 Bio	6 x 1 ml	08 083 819 001
Control E Level 3 Bio	6 x 1 ml	08 083 827 001

## Regulatory disclaimer

For use in quality control/manufacturing process only.

## Trademarks

CEDEX is a trademark of Roche.

## Published by

Roche Diagnostics GmbH  
Sandhofer Straße 116  
68305 Mannheim  
Germany

© 2020 Roche Diagnostics GmbH  
All rights reserved.