

NxtScript enzyme and master mix

Next generation high performance RT-PCR

NxtScript RT for high temperature transcription

NxtScript RT is a designed M-MLV reverse transcriptase with high thermostability. The enzyme can be used to generate first-strand cDNA from mRNA or total RNA for use in down-stream applications, such as two-step RT-PCR, RT-PCR for detection of viral and mRNA targets (e.g., cancer biomarkers), and generation of full-length cDNA libraries. Point mutations increase the thermostability of the enzyme and support higher cDNA yield of full-length transcripts than wild type M-MLV reverse transcriptase.

NxtScript DNA Master for stable multiplex performance

The NxtScript DNA Master is a sensitive and robust 5x conc. master mix. You can use it to run qPCR reactions and detect DNA, or use it in combination with NxtScript Reverse Transcriptase to detect RNA targets, such as viral pathogens. This mix is optimized for multiplex reactions and shows high room temperature stability.

- **Easily detect difficult DNA targets**



The high thermostability of NxtScript allows reactions up to 60°C to overcome RNA secondary structures (e.g., in GC-rich templates).

- **Save time and stay specific**

Take advantage of the aptamer hot start technology and run a fast PCR protocol.

- **Benefit from high quality standards and multiple formats**

Roche quality management is inserted in every production step for full documentation and traceability. Customized manufacturing and filling upon request.



Reverse transcription at different temperatures and subsequent qPCR show the high thermostability of NxtScript. Reference: data on file at Roche.

Ordering information

Product	Pack size	Catalog Number
NxtScript RT, 85 U/μl, 80 μl	80 μl	07 371 527 103
NxtScript Reverse Transcriptase, conc.	custom fill	07 051 166 103

Specification

Appearance	clear, colorless solution
Activity	750 – 2200 U/μl
Purity (HPLC)	≥90%
Nonspecific endonucleases (MWM III DNA)	not detectable with up to 75 U after 16 hours incubation at 37°C
Ribonucleases (MS2 RNA)	not detectable with up to 150 U after 1 hour incubation at 37°C
Stability	12 months at –15 to –25°C

Product	Pack size	Catalog Number
NxtScript DNA Master	5 ml	07 368 143 103

Specification

Appearance	clear, colorless solution
Function tested in qPCR (human reference cDNA G6PDH/β2M fragments)	corresponds
Stability	37 months at –15 to –25°C

Related products	Pack size	Catalog Number
NxtScript Reverse Transcriptase, high conc., 1MU	200 μl	09 917 225 001
NxtScript 2G RT, conc.	custom fill	09 085 220 103
AptaTaq DNA Polymerase, 5 U/μl	custom fill	05 457 882 103
AptaTaq DNA Polymerase, glycerol-free, 50 U/μl	custom fill	05 187 605 103
KAPA3G HotStart DNA Polymerase, glycerol-free, 30 U/μl	custom fill	08 918 651 103
KAPA3G HotStart Master	custom fill	09 084 711 103

© 2024
All rights reserved.

Published by

Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany

[custombiotech.roche.com](https://www.custombiotech.roche.com)

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America
mannheim.custombiotech@roche.com

United States
custombiotech.ussales@roche.com

Canada
custombiotech.can@roche.com

Asia Pacific
apac.custombiotech@roche.com

Regulatory Disclaimer
For further processing only.