

Carboxypeptidase B, recombinant in *Pichia pastoris*

New standards in secure industrial applications

Safer products, starting from the source

As industrial production of biological reagents and biopharmaceuticals expands, heightened regulatory requirements aim to ensure safer, cleaner products. Eliminating sources of animal material within your production process is a simple and effective way of complying with stringent regulations, but equally important is to maintain robust and optimized manufacturing. CustomBiotech has transformed the rat carboxypeptidase B gene into the expression host *Pichia pastoris* to create a recombinant enzyme that meets evolving regulatory demands and delivers high performance.

Optimal for sensitive manufacturing workflows

Carboxypeptidase B provides critical process safety for biopharma workflows, such as active insulin production and post-translational modification of antibodies.

- Animal-component-free production, reducing risk of contamination
- Manufactured under DIN EN ISO 13485 in accordance to good manufacturing practice (GMP) guidelines

Consistency that saves time and money

Each enzyme lot meets narrow quality and activity specifications so you can count on not only a quality supply, but also consistently high performance.

- Supports robust production processes
- Less downtime and corrective action required
- High reproducibility of your own product

A smart addition to standing operations

Carboxypeptidase B is easily integrated into standing manufacturing workflows, leveraging Roche expertise and support.

- Dedicated support team can assist in optimizing production steps and ensuring a supply that meets your specific needs
- Regulatory expertise available upon request to support product registration

Specification

Carboxypeptidase B from rat, recombinant in *Pichia pastoris*

Appearance	clear, colorless to slightly yellowish solution
Composition	Tris/HCl, 33 mmol/l; ZnCl ₂ , 0.1 mmol/l
pH value	7.5 – 8.5
Activity	≥ 400 U/ml (with Hippurylarginine)
Specific activity	≥ 210 U/mg (with Hippurylarginine)
Protein	2.8 ± 0.5 mg/ml
Purity (RP_HPLC)	≥ 85%
Bioburden	≤ 50 cfu/ml
Trypsin activity	≤ 0.005%
Stability	≥ 12 months at -15 to -25 °C

Ordering information

Product	Pack size	Catalog Number
Carboxypeptidase B from rat, recombinant in <i>Pichia pastoris</i>	30 KU	03 358 682 103
Related products	Pack size	Catalog Number
Trypsin, recombinant from <i>Pichia pastoris</i>	3.5 MU	03 358 658 103
Proteinase K, recombinant PCR Grade	850 ml	03 654 672 103

Regulatory Disclaimer

For further processing only.

© 2023
All rights reserved.

Published by

Roche Diagnostics GmbH
Sandhofer Str. 116
68305 Mannheim
Germany

[custombiotech.roche.com](https://www.custombiotech.roche.com)

Please contact your local CustomBiotech representative

Europe, Middle East, Africa, Latin America
mannheim.custombiotech@roche.com

United States
custombiotech.ussales@roche.com

Canada
custombiotech.can@roche.com

Asia Pacific
apac.custombiotech@roche.com